

[Food Hydrocolloids, 4, 323 (1990)]

The Decrease of Thaumatin's Sweetness Intensity Upon Interaction with Carrageenan.

SHIRO OHASHI, FUMIKO URA, MASANORI TAKEUCHI, HIROKI IIDA,
KAZUSHI SAKUE, TAKASHI OCHI, SHIGEO UKAI*, KOICHI HIRAMATSU

Carrageenan was added at various ratios to thaumatin, a sweet protein, and the interactions between thaumatin and λ -, κ - and ι -carrageenan were investigated from the following standpoints: pH, turbidity at 550nm, CD spectral change ($\Delta\epsilon$) and the decrease in sweetness intensity. Decrease in thaumatin sweetness intensity was observed with lower pH (<4) and higher carrageenan concentrations in the following order: $\lambda \approx \iota > \kappa$ -carrageenan. There was no correlation between sweetness intensity decrease and turbidity. With the increase in carrageenan concentration, CD spectral change ($\Delta\epsilon$) was distinctly observed at pH 3-4. Thaumatin's sweetness intensity decreased simultaneously.

[Food Hydrocolloids, 4, 379 (1990)]

Interaction of Thaumatin with Carrageenans. II. Effects of pH, Temperature and Competing Cations studied by Circular Dichroism.

SHIRO OHASHI, FUMIKO URA, MASANORI TAKEUCHI, HIROKI IIDA,
KAZUSHI SAKAUE, TAKASHI OCHI, SHIGEO UKAI*, KOICHI HIRAMATSU

Interaction of thaumatin with carrageenans was studied spectroscopically on the basis of circular dichroism (CD) at pH 3-7 and 20-60°C and examined from estimating the sweetness of the complex. The most noticeable reduction in sweetness intensity of thaumatin was pH 3-4 and less sweetness reduction occurred at pH 5-7. Thaumatin was completely dissociated from a thaumatin-carrageenan complex at pH 5-7 when a salt was added to the complex, but no complete dissociation occurred at pH 3-4. The results suggest that the complex formation at pH 5-7 is due to electrostatic bonding (ionic bonds), but that at pH 3-4 another bonding mechanism is also involved.

[Free Radical Biol. Med., 8, 25 (1990)]

A One-step Enzyme Immunoassay for Human Manganese Superoxide Dismutase with Monoclonal Antibodies.

TETSUO ADACHI*, KAZUYUKI HIRANO, KYOZO HAYASHI, YASUTOSHI MUTO,
FUMITAKA OKUNO

A one-step enzyme immunoassay for the determination of manganese superoxide dismutase in serum has been developed with two kinds of monoclonal antibodies. Proposed method had high sensitivity (assay range, 0.4~200 ng/ml), good recovery (recovery percentage, 102.9~106.2%) and reproducibilities (intraassay, C.V.=1.87~3.66%; interassay, C.V.=3.03~10.4%). From these results, it is possible to apply this method to routine clinical analysis and biochemical research with various purposes.